



TRANSPARENT MICRO-ADHERENT WOUND CONTACT LAYER WITH SAFETAC


- Atraumatic to the wound and surrounding skin on removal
- Can remain in place for up to 14 days which promotes cost-effective wound healing

The dressing

- Medical grade silicone bound to an elastic polyamide net
- Non-absorbent
- Conforms to body contours
- Transparent

Release film

- Protects the Safetac layer
- Allows sterile application



Safetac® layer (soft silicone)

- Non-adherent to moist wound bed
- Adheres gently to surrounding intact skin
- Minimizes trauma and pain at dressing change



Safetac® technology

Safetac® technology. Less trauma. Less pain.™

Safetac is a patented soft silicone adhesive technology. Dressings with Safetac technology are atraumatic both during wear and upon removal. These dressings minimize trauma to the wound and the surrounding skin, which minimizes pain to the patient. They also prevent maceration by sealing around the margin of the wound to protect peri-wound skin.



Skin stripping occurs with traditional adhesive




No skin stripping occurs with Safetac


How Mepitel works

Mepitel can be left in place for up to 14 days depending on the condition of the wound. The porous structure of Mepitel allows exudate to pass into an outer absorbent dressing. Mepitel with Safetac soft silicone adhesive technology prevents the outer dressing from sticking to the wound and therefore minimizes trauma and pain associated with dressing changes. This results in less discomfort for the patient and which in turn ensures undisturbed wound healing.


How to use Mepitel



Gently cleanse the wound; dry surrounding skin. Remove the release film.



Apply Mepitel directly on to the wound allowing it to overlap onto the surrounding skin by 2cm. Remove the remaining backing.



Apply a cover dressing over Mepitel and fixate in place.

Precautions

When used on bleeding wounds or wounds with high viscosity exudate, Mepitel should be covered with a moist absorbent dressing pad.

When used for fixation of skin grafts, on burns with blisters or for skin tears, Mepitel should be left in place for 5 to 7 days.

Mepitel



Areas of use

Mepitel is designed for a wide range of wounds such as skin tears, skin abrasions, surgical incisions, second degree burns, blistering, lacerations, partial and full thickness grafts, diabetic ulcers, venous and arterial ulcers, pressure ulcers, digit injuries and Epidermolysis Bullosa (EB). Mepitel can also be used on wounds associated with Pediatrics where management of pain associated with dressing change is a primary consideration.

Benefits of Mepitel

- Does not adhere to a moist wound bed but adheres gently to the dry peri-wound area.
- Minimizes trauma and pain during dressing changes.
- Can remain in place for up to 14 days depending on wound condition.
- Seals around the wound margins to reduce the risk of maceration.
- Conforms well to body contours.
- Topical preparations can be applied as required, leaving Mepitel in place.
- Maintains an optimal moisture balance when used with appropriate cover.

Mepitel Assortment (Sterile packed)

Art. No	Size cm	Pcs/box	Pcs/case
290510	5 x 7.5	10	50
290710	7.5 x 10	10	40
291010	10 x 18	10	70
292005	20 x 30	5	30