

# Stays on


Uniquely conforms



## 360° Flex technology – stays on and uniquely conforms

- Improves comfort<sup>1,2</sup>
- Keeps the dressing in place<sup>1,2</sup>

## Y-cuts expand to allow movement



## Safetac® wound contact layer

- Less painful dressing changes<sup>6</sup> and reduced risk of maceration<sup>7-8</sup>



## Highly breathable backing film

- Exudate Progress Monitor lets you track and record exudate progress<sup>3</sup>
- Backing film provide a barrier to microorganisms larger than 25nm<sup>4</sup>
- Water resistant so your patients can shower<sup>5</sup>



## Spreading layer

- Distributes exudate evenly over a wide surface area to maximize exudate management<sup>9</sup>



## Foam layer

- Absorbs exudate and transports it to the spreading layer<sup>9</sup>
- The 4-layer construction works to absorb and channel exudate away from the wound bed<sup>9</sup>

# Mepilex® Border Flex Lite

Our next generation of flexible dressings



## Designed to conform and stay on

- Innovative 360° Flex Technology allows Mepilex® Border Flex Lite to move in every direction, reducing skin stress, increasing comfort and keeping it in place<sup>1,2</sup>
- Up to 77% more conformable than Mepilex Border Lite<sup>1,2</sup>

## Smart exudate management

- The unique 4-layer dressing absorbs and channels exudate away from the wound bed<sup>3</sup> and allows you to track exudate progress without disturbing the wound<sup>9</sup>

## Safetac® technology

- Our proven Safetac Technology means less pain and trauma to the skin during dressing changes<sup>6</sup> and less maceration compared to traditional (adhesive) dressings.<sup>7,8</sup>

## Areas of use

For the management of non-to-low exuding wounds, including:

- Leg and foot ulcers
- Pressure injuries
- Surgical and traumatic wounds such as abrasions, blisters and skin tears
- May be used for protection of compromised and/or fragile skin

## Easy application




The largest size of Mepilex Border Flex Lite is equipped with our new proprietary three-part release line.

\*Available in 15x15 cm

**Note:** If you see signs of clinical infection e.g. fever or if the wound or surrounding skin becomes red, warm or swollen, consult a health care professional for appropriate treatment.

## Ordering information (sterile packed)

| Product Code | Product size | Pcs/Box | Pcs/Case |
|--------------|--------------|---------|----------|
| 581011       | 4 x 5 cm     | 10      | 70       |
| 581100       | 5 x 12.5 cm  | 5       | 65       |
| 581200       | 7.5 x 7.5 cm | 5       | 70       |
| 581300       | 10 x 10 cm   | 5       | 50       |
| 581500       | 15 x 15 cm   | 5       | 50       |



**References:** **1.** Mölnlycke Health Care. PD-570458. Data on file. **2.** Mölnlycke Health Care. Comparison between MBFL and AGBL. Technical Report 002M. Data on file. **3.** Mölnlycke Health Care. PD-570599. Data on file. **4.** Mölnlycke Health Care. PD-580421. Data on file. **5.** Mölnlycke Health Care. PD-580483. Data on file. **6.** White, R.A. A multinational assessment of pain when removing dressings. Wounds UK 2008. 4(1): 14-22. **7.** Woo, K. & Bergström, C. An Open Randomized Cross-over Investigation Assessing Perceived Pain at Dressing Change Comparing A Soft Silicone Dressing, Mepilex® Border, with an Adhesive hydrocellular polyurethane dressing, Allevyn® Adhesive in Patients with chronic ulcer [Clinical Investigation Report (MXB408, Part A, Mölnlycke Health Care), 2007]. **8.** Meaume, S et al. A Study to Compare a New Self-Adherent Soft Silicone Dressing with a Self-Adherent Polymer Dressing in Stage II Pressure Ulcers. Ostomy Wound Manag. 2003; 49(9): 44-51. **9.** Mölnlycke Health Care. PD-570117. Data on file.